

# Moving Cryosurgery into the 21<sup>st</sup> Century!



## cryopen®


The CryoPen Cryosurgical System provides **Simple, Safe** and **Effective** cryosurgical treatment for common skin lesions without dangerous cryogenic gases or liquids. The CryoPen enables doctors to dramatically increase revenue while eliminating the need for expensive spray canisters typically required for cryogenic procedures .

- **Simple.** Non-technique dependent procedure for treatment of common skin lesions
- **Safe.** No more handling of dangerous cryogenic gases or liquids such as liquid nitrogen
- **Effective.** Pen-point precision with consistent freeze temperatures of -105 Celsius

**Significant Return on Investment.**  
Most practices recover their entire investment in less than 6 months!

CryoPen.com | 1-888-246-3928

The CryoPen Cryosurgical System was designed to replace outdated cryosurgical instruments by using the newest cooling technology to eliminate the need for cryogenic gases or liquids.



## The CryoPen

The CryoPen uses no cryogenic gases or liquids. It offers the most advanced state-of-the-art cooling technology in the field of cryosurgery. The CryoPen Cryosurgical System is designed to be the most effective and simplest cryosurgical instrument on the market. It is less costly than other current treatments, and is less painful than traditional methods. The CryoPen also offers accurate pen-point precision and a consistent freeze temperature, which can be monitored using the CryoPen Temperature Indicator. With its simple, non-technique dependent procedure, the CryoPen makes cryosurgery safe, simple, and effective.

## Common Lesions Treated

- warts/plantar warts
- actinic keratoses
- lentigo
- keloids
- dermatafibroma
- skin tags
- seborrheic keratoses
- non-recurrent basal cell
- molluscum

### Patient Benefits

- No exposure to dangerous cryogenic gases or liquids
- Immediate treatment without referral
- No anesthetic needed
- Less painful than other methods currently available

### Physician Advantage

- Efficacy and consistent patient outcomes
- Incredibly easy to use
- Increased patient convenience, safety, and satisfaction

### Practice Advantage

- Significant revenue generator and profit center
- No hazardous gases, liquids, or OSHA requirements
- Keeps patient care and reimbursements in-house

### Breakthrough Technology

**Temperature Indicator:**  
Allows for verification of consistent, optimal treatment temperature to maximize efficacy and minimize variations between procedures.

### Chilling Wells:

With state-of-the-art cooling technology, each CryoPen reaches the optimal surgical temperature of -105°C without cryogenic gases or liquids.

### CryoPen:

The CryoPen cryosurgical system provides four handheld CryoPens.

### Re-Usable Tips:

CryoPen re-usable tips ensure precise treatment areas and are available in sizes of 3mm, 5mm, 7mm, and 10mm.



cryopen.

**Contact us for more information**

1-888-246-3928 • [www.cryopen.com](http://www.cryopen.com) • [sales@cryopen.com](mailto:sales@cryopen.com)

All Contents Copyright CryoPen® Inc. 2008-2011 All Rights Reserved.